


Monoclonal Antibody to Ferret IFN- γ

PF647P CONJUGATED

Antibody:	MTF14
Product code:	3112-72-100T
Size:	100 tests
Immunogen:	Recombinant ferret IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant ferret IFN- γ
Contents:	Ready-to-use solution of conjugated antibody in sterile filtered (0.2 μ m) PBS with 0.2% BSA and 0.09% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugation:	The MTF14 antibody has been conjugated to PF647P.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	Detection of ferret IFN- γ producing cells by flow cytometry. 5 μ l is recommended for staining of 1 million cells in a total volume of 50 μ l.



Detection of IFN- γ by flow cytometry in viable ferret spleen cells. Cells were stimulated for 16 hours in the presence of PMA/ ionomycin and Brefeldin A. Cells were then fixed and permeabilized using 4% paraformaldehyde and saponin, and subsequently stained with MTF14-PF647P (right panel). The left panel depicts cells stained with a matched isotype control antibody. The dot plots derive from gated events of typical lymphocyte characteristics in forward and side light scatter. Flow cytometry was performed on a BD FACSVerse system.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages there from.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2022-11-14

www.mabtech.com