

Monoclonal Antibody to Bovine Interleukin-8 (CXCL8)

BIOTINYLATED

Antibody:	26E5
Product code:	3114-6-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-8
Isotype:	Mouse IgG2b
Specificity:	Native and recombinant bovine IL-8. Native and recombinant human IL-8. Cross-reacts with IL-8 from cynomolgus macaques and IL-8 from dog. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-8 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-8 producing cells using ELISpot. 26E5 is recommended as detection mAb: in ELISA in combination with capture mAb MT8H6 (product code 3114-3) and in ELISpot in combination with coating mAbs MT8H6/8F19 (product code 3114-1).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com