

# Monoclonal Antibody to Bovine Interferon- $\gamma$

BIOTINYLATED

<b>Antibody:</b>	MT307
<b>Product code:</b>	3119-6-250
<b>Quantity:</b>	250 $\mu$ g
<b>Immunogen:</b>	Recombinant bovine IFN- $\gamma$
<b>Isotype:</b>	Mouse IgG2a
<b>Specificity:</b>	Native and recombinant bovine IFN- $\gamma$ . Cross-reacts with native IFN- $\gamma$ from sheep and goat.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IFN- $\gamma$ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- $\gamma$ producing cells using ELISpot. MT307-biotin is recommended as detection mAb in combination with coating mAb MT17.1 (product code 3119-3). MT307 is also recommended for use in flow cytometry.
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on, 2022-11-18

[www.mabtech.com](http://www.mabtech.com)