

Monoclonal Antibody to Porcine Interferon- γ

BIOTINYLATED

Antibody:	P2C11
Product code:	3130-6-1000
Quantity:	1000 μ g
Immunogen:	Recombinant porcine IFN- γ
Isotype:	Mouse IgG2a
Specificity:	Native and recombinant porcine IFN- γ
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C.
Applications:	For quantification of porcine IFN- γ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- γ producing cells using ELISpot. P2C11 is recommended as detection mAb in combination with coating mAb pIFN γ -I (product code 3130-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com