

Monoclonal Antibody to Porcine TNF- α

BIOTINYLATED

Antibody:	MT859
Product code:	3135-6-1000
Quantity:	1000 μ g
Immunogen:	Recombinant porcine TNF- α
Isotype:	Mouse IgG2a
Specificity:	Native and recombinant porcine TNF- α
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of porcine TNF- α in solution e.g. cell culture supernatants using ELISA and for enumeration of porcine TNF- α producing cells using ELISpot. MT859 is recommended as detection mAb in ELISA in combination with coating mAb MT1P14 (product code 3135-1) and in ELISpot in combination with coating mAb MT32E9 (product code 3135-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com