

# Monoclonal antibody to Bovine IgG

<b>Antibody:</b>	MT134
<b>Product code:</b>	3150-3-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Bovine IgG
<b>Isotype:</b>	Rat IgG1
<b>Specificity:</b>	Bovine IgG Cross-reacts with IgG from human (86% cross-reactivity), horse (6% cross-reactivity) and sheep (1% cross-reactivity).
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For enumeration of bovine IgG producing cells using ELISpot. MT134 is recommended as coating mAb in combination with detection mAb MT391 (product code 3150-6). MT134 can also be used as capture mAb in ELISA in combination with detection mAb MT391.
<b>Protocols:</b>	Please contact Mabtech for ELISpot protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com