

Monoclonal Antibody to Porcine IgG

BIOTINYLATED

Antibody:	MT424
Product code:	3151-6-1000
Quantity:	1000 µg
Immunogen:	Porcine IgG
Isotype:	Mouse IgG2a
Specificity:	Porcine IgG No cross-reactivity (<0.1%) with IgG from cow, horse, sheep, goat, human and rat.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For enumeration of porcine IgG producing cells using ELISpot. MT424 is recommended as detection mAb in combination with coating mAb MT421 (product code 3151-3). MT424 can also be used as detection mAb in ELISA in combination with coating mAb MT421.
Protocols:	Please contact Mabtech for ELISpot protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com