

Monoclonal antibody to Rabbit IgG

Antibody:	MT5131
Product code:	3152-3-1000
Quantity:	1000 µg
Immunogen:	Rabbit IgG
Isotype:	Mouse IgG1
Specificity:	Rabbit IgG No cross-reactivity (<0.1%) with rabbit IgA, IgM, or with IgG from human, monkey, mouse, guinea pig, cow, horse, pig, dog, and cat.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IgG in solution e.g. serum/plasma in ELISA and for enumeration of IgG producing cells using ELISpot. MT5131 is recommended as coating mAb in combination with detection mAb MT536 (product code 3152-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com