

Monoclonal Antibody to Rabbit IgG

BIOTINYLATED

Antibody:	MT536
Product code:	3152-6-1000
Quantity:	1000 µg
Immunogen:	Rabbit IgG
Isotype:	Mouse IgG1
Specificity:	Rabbit IgG No cross-reactivity (<0.1%) with rabbit IgA, IgM, or with IgG from human, monkey, mouse, guinea pig, cow, horse, pig, dog, and cat.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IgG in solution e.g. serum/plasma in ELISA and for enumeration of IgG producing cells using ELISpot. MT536 is recommended as detection mAb in combination with coating mAb MT5131 (product code 3152-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com