

# Monoclonal antibody to Equine IgG

<b>Antibody:</b>	MT329A4
<b>Product code:</b>	3153-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Equine IgG
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Equine IgG No cross-reactivity (<0.01%) with equine IgA, IgM, or with IgG from cow, sheep, goat, pig, dog, cat, ferret, guinea pig, hamster, rabbit, rat, mouse, macaques, and human.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of equine IgG in solution e.g. serum/plasma in ELISA. MT329A4 is recommended as coating mAb in combination with detection mAb MT310A61 (product code 3153-6).
<b>Protocols:</b>	Please contact Mabtech protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com