

Goat anti-mouse IgG, HRP

Product code: 3310-14-1000

Antibody: Antiserum is raised in goat using purified antigen as immunogen. The antibody is specific for mouse IgG H-chains (Fc_γ) and was purified using immunoaffinity chromatography. No reactivity against human, horse and bovine serum proteins or against mouse IgM, IgA or non-immunoglobulin serum proteins. Conjugate cross-reacts with rat immunoglobulins and may cross-react with immunoglobulins from other species.

Enzyme: Horseradish peroxidase

Conjugation: The antibody is conjugated to highly purified enzyme from horseradish. The conjugate has been thoroughly depleted of free ligand and free enzyme.

Quantity: 1 ml

Applications: ELISA, Western blot

Working dilution: Supplied in storage buffer with 0.002% Kathon CG. Dilute 1/1000 as a general guideline. Dilution factor may vary with differing assay conditions.

Storage: Shipped at ambient temperature. Store at 4-8 °C upon receipt.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com