

Monoclonal Antibody to Mouse Interferon- α

BIOTINYLATED

Antibody:	MT9L/14A/113
Product code:	3326-14-1000
Quantity:	1000 μ g
Immunogen:	Recombinant mouse IFN- α
Isotype:	Rat IgG1
Specificity:	Native and recombinant mouse IFN- α . MT9L/14A/113 in combination with antibodies MT24A/44A/104 detect subtypes 1, 2, 4, 5, 6T, 7/10 (129/Sv), 8/6, 9, 11, 12, 13, 14, 15 (A), and B.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IFN- α in solution e.g. cell culture supernatants using ELISA. MT9L/14A/113 are recommended as detection mAbs in combination with coating mAbs MT24A/44A/104 (product code 3326-1).
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com