

# Monoclonal Antibody to Human IP-10

BIOTINYLATED

<b>Antibody:</b>	MTP13
<b>Product code:</b>	3401-6-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Recombinant human IP-10
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant human IP-10 (CXCL10). Cross-reacts with IP-10 from rhesus and cynomolgus macaques. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For enumeration of IP-10 producing cells using ELISpot. MTP13 is recommended as detection mAb in combination with coating mAb MTP11 (product code 3401-3).
<b>Protocols:</b>	Please contact Mabtech for protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)

Updated on, 2022-11-14

[www.mabtech.com](http://www.mabtech.com)