


Monoclonal Antibody to Human IL-4

PF488P CONJUGATED

- Antibody:** IL4-3
- Product code:** 3410-71-100T
- Size:** 100 tests
- Immunogen:** Recombinant human IL-4
- Isotype:** Rat IgG1
- Specificity:** Human IL-4. Cross-reacts with IL-4 from cynomolgus macaque. Please visit www.mabtech.com for information on reactivity with other monkey species.
- Contents:** Ready-to-use solution of conjugated antibody in sterile filtered (0.2µm) PBS with 0.2% BSA and 0.09% sodium azide.
- Purification:** Purified from *in vitro* cultures by protein G affinity chromatography.
- Conjugation:** The antibody has been conjugated to PF488P.
- Storage:** Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
- Applications:** Detection of IL-4 producing cells by flow cytometry. 5 µl is recommended for staining of 1 million cells in a total volume of 50 µl.



Detection of IL-4 by flow cytometry in viable human PBMC. Cells were stimulated for 16 hours in the presence of PMA/ionomycin and Brefeldin A. Cells were stained with anti-CD3-BV421, fixed and permeabilized using 4% paraformaldehyde and saponin, and further stained with PF488P-conjugated mAb IL4-3 (right panel). Non-stimulated cells shown as reference (left panel). The dot plot derives from gated events of typical lymphocyte characteristics in forward and side light scatter. Flow cytometry was performed on a BD FACSVerse system.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages there from.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2022-11-14

www.mabtech.com