

Monoclonal Antibody to Human Interferon- γ

BIOTINYLATED

Antibody:	MT111W
Product code:	3420-2-250
Quantity:	250 μ g
Immunogen:	Recombinant human IFN- γ
Isotype:	Rat IgG2a
Specificity:	Native and recombinant human IFN- γ . Native and recombinant IFN- γ from cynomolgus macaque. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	Western blot. The antibody is also recommended for neutralization of human IFN- γ bioactivity. For neutralization assays, mAb MT111W is available unconjugated and free of preservatives (product code 3420-1N-500).
Protocols:	Please contact Mabtech for more information.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com