

Monoclonal Antibody to Human IFN- γ

Antibody:	1-D1K
Product code:	3420-7
Quantity:	100 tests
Immunogen:	Recombinant human IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IFN- γ
Contents:	Ready-to-use solution of conjugated antibody in sterile filtered (0.2 μ m) PBS with 0.2% BSA and 0.09% sodium azide. 5 μ l is recommended for staining of 1 million cells in a total volume of 50 μ l.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugation:	The 1-D1K has been conjugated to fluorescein isothiocyanate (FITC).
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	Detection of IFN- γ producing cells by immunofluorescent staining of cells or tissue sections.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2013-11-19



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com