

Monoclonal Antibody to Human Interleukin-3

BIOTINYLATED

Antibody:	IL3B
Product code:	3435-6-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-3
Isotype:	Rat IgG1
Specificity:	Native and recombinant human IL-3
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-3 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-3 producing cells using ELISpot. IL3B is recommended as detection mAb in combination with coating mAb IL3A (product code 3445-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com