

# Monoclonal antibody to Human Interleukin-2

<b>Antibody:</b>	MT2A91
<b>Product code:</b>	3445-1-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Recombinant human IL-2
<b>Isotype:</b>	Mouse IgG2b
<b>Specificity:</b>	Native and recombinant human IL-2. Native and recombinant IL-2 from cynomolgus macaque. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	Western blot
<b>Protocols:</b>	Please contact Mabtech for more information.

# MABTECH

Note: for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)