

Monoclonal antibodies to Human Interleukin-2

Antibody:	MT2A91/2C95 The reagent consists of a combination of two different monoclonal antibodies.
Product code:	3445-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-2
Isotype:	Mouse IgG2b/Mouse IgG1
Specificity:	Native and recombinant human IL-2. Native and recombinant IL-2 from cynomolgus macaque. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-2 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-2 producing cells using ELISpot. MT2A91/2C95 are recommended as coating mAbs in combination with detection mAb MT8G10 (product code 3445-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com