

Monoclonal antibodies to Human Interleukin-27

Antibody:	MT27/361
Product code:	3458-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-27
Isotype:	Mouse IgG2b/IgG1
Specificity:	Native and recombinant human IL-27. MT27 binds to IL-27 and MT361 binds to p28, which is part of the IL-27 heterodimer.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-27 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-27 producing cells using ELISpot. MT27/361 are recommended as coating mAbs in combination with detection mAb MT140 (product code 3458-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com