

Monoclonal Antibody to Human EBI3

BIOTINYLATED

Antibody:	MT140
Product code:	3458-6-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-27
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human EBI3. MT140 binds to the EBI3 subunit of IL-27, IL-35 and IL-39.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-27 and EBI3 in solution e.g., cell culture supernatants using ELISA and for enumeration of IL-27 and EBI3 producing cells using ELISpot. MT140 is recommended as detection mAb for analysis of: <ul style="list-style-type: none">- IL-27 in combination with coating mAbs MT27/361 (product code 3458-3)- EBI3 in combination with coating mAb MT23A14 (product code 3459-3)
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-18

www.mabtech.com