

Monoclonal Antibody to Human Interleukin-6

BIOTINYLATED

Antibody:	39C3
Product code:	3460-6-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-6
Isotype:	Rat IgG2a
Specificity:	Native and recombinant human IL-6. Cross-reacts with native IL-6 from rhesus and cynomolgus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-6 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-6 producing cells using ELISpot. 39C3 is recommended as detection mAb in combination with coating mAb 13A5 (product code 3460-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com