

Monoclonal antibody to Monkey/Human Interleukin-13

Antibody:	IL13-I
Product code:	3470M-3-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-13
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IL-13. Cross-reacts with native IL-13 from rhesus and cynomolgus macaques and common marmosets. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-13 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-13 producing cells using ELISpot. IL13-I is recommended as coating mAb in combination with detection mAb IL13-3 (product code 3470-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com