

# Monoclonal antibody to Human Interleukin-13

<b>Antibody:</b>	25K2
<b>Product code:</b>	3471-5-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Recombinant human IL-13
<b>Isotype:</b>	Rat IgG1
<b>Specificity:</b>	Native and recombinant human IL-13. Cross-reacts with native IL-13 from rhesus and cynomolgus macaques and common marmosets. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	See vial. The antibody is supplied in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IL-13 in solution e.g. cell culture supernatants using ELISA. 25K2 is recommended as coating mAb in combination with detection mAb MT1318 (product code 3471-4).
<b>Protocols:</b>	Please contact Mabtech for ELISA protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)