

# Monoclonal Antibody to Human GM-CSF

BIOTINYLATED

<b>Antibody:</b>	23B6
<b>Product code:</b>	3480-6-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Recombinant human GM-CSF
<b>Isotype:</b>	Rat IgG2a
<b>Specificity:</b>	Native and recombinant human GM-CSF. Cross-reacts with native GM-CSF from rhesus and cynomolgus macaques. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of GM-CSF in solution e.g. cell culture supernatants using ELISA and for enumeration of GM-CSF producing cells using ELISpot. 23B6 is recommended as detection mAb in combination with coating mAb 21C11 (product code 3480-3).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note: for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)

Updated on, 2022-11-14

[www.mabtech.com](http://www.mabtech.com)