

Monoclonal Antibody to Human Granzyme B

BIOTINYLATED

Antibody:	MT28
Product code:	3486-4-1000
Quantity:	1000 µg
Immunogen:	Recombinant human Granzyme B
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human Granzyme B
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of Granzyme B in solution e.g. cell culture supernatants using ELISA and for enumeration of Granzyme B producing cells using FluoroSpot. MT28 is recommended as detection mAb in combination with coating mAb MT34B6 (product code 3486-1).
Protocols:	Please contact Mabtech for ELISA and FluoroSpot protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com