

Monoclonal antibody to Monkey Granzyme B

| | |
|-----------------------|---|
| Antibody: | MTAG8 |
| Product code: | 3486M-3-1000 |
| Quantity: | 1000 µg |
| Immunogen: | Recombinant rhesus macaque granzyme B |
| Isotype: | Mouse IgG1 |
| Specificity: | Native and recombinant granzyme B from rhesus macaque and cynomolgus macaque. Please visit www.mabtech.com for information on reactivity with other monkey species. |
| Concentration: | Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide. |
| Purification: | Purified from <i>in vitro</i> cultures by protein G affinity chromatography. |
| Storage: | Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use. |
| Applications: | For quantification of granzyme B in solution e.g. cell culture supernatants using ELISA and for enumeration of granzyme B producing cells using ELISpot. MTAG8 is recommended as coating mAb in combination with detection mAb MT8213 (product code 3486M-6). |
| Protocols: | Please contact Mabtech for protocols. |

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com