

Monoclonal antibody to Human CCL2

Antibody:	MTBH5
Product code:	3491-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant human CCL2 (MCP-1)
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human CCL2 (MCP-1). Native CCL2 (MCP-1) from rhesus and cynomolgus macaque. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of CCL2 (MCP-1) in solution e.g. cell culture supernatants using ELISA. MTBH5 is recommended as coating mAb in combination with detection mAb MTBH9-biotin (product code 3491-6).
Protocols:	Please contact Mabtech for protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com