

Monoclonal antibody to Human CCL4

BIOTINYLATED

Antibody:	MTBN3
Product code:	3494-6-1000
Quantity:	1000 µg
Immunogen:	Recombinant human CCL4 (MIP-1β)
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human CCL4 (MIP-1β). Native CCL4 (MIP-1β) from rhesus and cynomolgus macaque. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of CCL4 (MIP-1β) in solution e.g. cell culture supernatants using ELISA. MTBN3 is recommended as detection mAb in combination with coating mAb MTBN13 (product code 3494-3).
Protocols:	Please contact Mabtech for protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2024-05-08

www.mabtech.com