

# Monoclonal antibody to Monkey/Human Interleukin-17A

<b>Antibody:</b>	MT241
<b>Product code:</b>	3520M-3-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Recombinant human IL-17A
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant human IL-17A. Native IL-17A from rhesus and cynomolgus macaques. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IL-17A in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-17A producing cells using ELISpot. MT241 is recommended as coating mAb in combination with detection mAb MT504 (product code 3520-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)