

Monoclonal Antibody to Human Interleukin-17F

BIOTINYLATED

Antibody:	MTF411
Product code:	3523-4-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-17F
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IL-17F and IL-17A/F. The antibody binds to the F subunit of the IL-17A/F heterodimer.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-17A/F in solution e.g. cell culture supernatants using ELISA. MTF411 is recommended as detection mAb in combination with coating mAb MT504 (product code 3520-5).
Protocols:	Please contact Mabtech for protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com