

Human Interleukin-31 Monoclonal antibody

FOR NEUTRALIZATION

General

Antibody:	MT313
Product code:	3530-0N-500
Quantity:	500 µg
Immunogen:	Recombinant human IL-31
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IL-31
Concentration:	1 mg/ml
Buffer:	Supplied in sterile filtered (0.2 µm) PBS. Contains no preservatives. The endotoxin content in the product is below 10 EU/mg.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	Neutralization of human IL-31 bioactivity. It is recommended to establish the optimal antibody concentration for the assay system used.


Neutralization

Neutralization assay: Neutralization of human IL-31 induced IL-6 production (see graph).

To determine the neutralizing capacity and ND_{50} , mixtures of human recombinant IL-31 (50 ng/ml) and various concentrations of mAb MT313 were added to human aortic endothelial cells that had been pre-stimulated with recombinant human IFN- γ (10 ng/ml). Following overnight incubation, cell supernatants were quantified in IL-6 ELISA.

Neutralization dose 50 (ND50): MT313 neutralized the effect of human IL-31 at a ND_{50} value of 0.1-0.6 $\mu\text{g/ml}$ when 50 ng/ml of human IL-31 was present.

ND_{50} is defined as the concentration of antibody that results in 50% inhibition of the biological response under the experimental conditions specified.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.



The products are for research use only.
MABTECH shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages there from.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com