

Monoclonal Antibody to Human Interleukin-8 (CXCL8)

BIOTINYLATED

Antibody:	26E5
Product code:	3560-14-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-8
Isotype:	Mouse IgG2b
Specificity:	Native and recombinant human IL-8. Cross-reacts with IL-8 from cynomolgus macaque, IL-8 from cow and IL-8 from dog. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For enumeration of human and monkey IL-8 producing cells using ELISpot. For quantification of monkey IL-8 in solution e.g. cell culture supernatants using ELISA. For IL-8 ELISpot, 26E5 is recommended as detection mAb in combination with coating mAbs MT8H6/8F19 (product code 3560-1). For monkey IL-8 ELISA, 26E5 is recommended as detection mAb in combination with coating mAb MT8H6 (product code 3560-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-18

www.mabtech.com