

# Monoclonal Antibody to Human CD4

<b>Antibody:</b>	CD4-1
<b>Product code:</b>	3606-3
<b>Quantity:</b>	100 tests
<b>Immunogen:</b>	Stimulated human leukocytes
<b>Isotype:</b>	Mouse IgG2a
<b>Specificity:</b>	Recognizes human CD4
<b>Contents:</b>	Ready-to-use solution of conjugated antibody in sterile filtered (0.2µm) PBS with 0.2% BSA and 0.09% sodium azide. 5 µl is recommended for staining of 1 million cells in a total volume of 50 µl.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Conjugation:</b>	The CD4-1 antibody has been conjugated to rhodamine.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
<b>Applications:</b>	Immunocytochemistry; Flow cytometry

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.