

Monoclonal Antibody to Human Apolipoprotein B

Antibody:	LDL 11
Product code:	3715-7
Quantity:	250 µg
Immunogen:	Human LDL
Isotype:	Mouse IgG1
Specificity:	Native human LDL
Contents:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugation:	The LDL 11 has been conjugated to fluorescein isothiocyanate (FITC).
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing
Applications:	For intracellular staining of LDL using immunofluorescence. The simultaneous addition of fixative and a non-ionic detergent is recommended.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.