

Monoclonal Antibody to Human Amyloid- β

BIOTINYLATED

Antibody:	<i>bm</i> -AbetaN
Product code:	3740-6-1000
Quantity:	1000 μ g
Immunogen:	Human Amyloid- β
Isotype:	Mouse IgG1
Specificity:	Synthetic and natural human Amyloid- β . Detects different forms of the peptide (mono- and oligomers, protofibrils and fibrils). Binds to residues 3-8.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For detection of Amyloid- β using ELISA, immunohistochemistry and Western blot. Other applications include immunoprecipitation and affinity purification of Amyloid- β . For questions and further information please contact Mabtech.
Protocols:	<p><i>Recommended dilutions:</i></p> <p>The following dilutions should serve as a general guideline. Further optimization may be required depending on the assay conditions.</p> <p>Western blot: 1:5000</p> <p>Immunohistochemistry: 1:5000. Tested on PFA-fixed, paraffin-embedded and PFA-fixed, fresh frozen sections. Heat-induced antigen retrieval and Formic acid treatment are recommended.</p> <p>ELISA: 1:5000 (when used as detecting antibody in indirect ELISA)</p>

The antibody was developed as part of a collaboration between Mabtech AB and BioArctic Neuroscience AB

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com