

# Monoclonal antibody to Human IgE

|                       |   |
|-----------------------|---|
| <b>Antibody:</b>      | 107   |
| <b>Product code:</b>  | 3810-3-1000   |
| <b>Quantity:</b>      | 1000 µg   |
| <b>Immunogen:</b>     | Human IgE   |
| <b>Isotype:</b>       | Mouse IgG1  |
| <b>Specificity:</b>   | Human IgE   |
| <b>Concentration:</b> | Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.   |
| <b>Purification:</b>  | Purified from <i>in vitro</i> cultures by protein G affinity chromatography.  |
| <b>Storage:</b>       | Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.  |
| <b>Applications:</b>  | For quantification of IgE in serum, plasma and cell culture supernatants using ELISA and for enumeration of IgE producing cells using ELISpot. 107 is recommended as coating mAb in combination with detection mAb 182 (product code 3810-6) for ELISA and with detection mAbs 107/182/101 (product code 3810-8) for ELISpot. |
| <b>Protocols:</b>     | Please contact Mabtech for ELISA protocols.   |

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com