

Goat anti-mouse IgG

BIOTINYLATED

Antibody:	Antiserum is raised in goat using purified antigen as immunogen.
Product code:	3825-6-250
Quantity:	250 µg
Immunogen:	Mouse IgG
Isotype:	NA
Specificity:	The antibody is specific for mouse IgG H-chains. No cross-reactivity with mouse IgA, IgM, human Ig or non-immunoglobulin serum proteins. Cross-reacts with rat immunoglobulins and may cross-react with immunoglobulins from other species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Immunoaffinity chromatography
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	Enzyme immunoassays, Western blot techniques
Protocols:	Please contact Mabtech for ELISpot protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-18

www.mabtech.com