

Monoclonal Antibody to Mouse IgG3

BIOTINYLATED

Antibody:	MG3
Product code:	3826-6-1000
Quantity:	1000 µg
Immunogen:	Mouse IgG3
Isotype:	Rat IgG1
Specificity:	Mouse IgG3. No cross-reactivity (<0.1%) with mouse IgG1, IgG2a, IgG2b, IgG2c (IgG2a analogue in e.g. C57BL/6 and NOD mice), IgA, IgE, IgM or with IgG from rat, goat, rabbit, cow and human.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For detection of mouse IgG3 in e.g. antigen-specific ELISA.
Protocols:	Please contact Mabtech for more information.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com