

Monoclonal antibody to Human IgA1

Antibody:	MT50
Product code:	3861-3-1000
Quantity:	1000 µg
Immunogen:	Human IgA
Isotype:	Rat IgG1
Specificity:	Human IgA1 derived from serum and colostrum. No cross-reactivity (< 0.1%) with human IgA2, IgG1, IgG2, IgG3, IgG4, IgE or IgM.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For detection of IgA1 in e.g. antigen-specific ELISA.
Protocols:	Please contact Mabtech for more information.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2022-11-17

www.mabtech.com