

# Monoclonal Antibody to Mouse IgA

BIOTINYLATED

<b>Antibody:</b>	MT39A
<b>Product code:</b>	3865-6-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Mouse immunoglobulins
<b>Isotype:</b>	Rat IgG1
<b>Specificity:</b>	Mouse IgA
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For enumeration of IgA (total and antigen-specific) producing cells using ELISpot and for quantification of antigen-specific IgA using e.g. indirect ELISA. In ELISpot, MT39A is recommended as detection mAb in combination with coating mAb MT45A (product code 3865-3).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note: for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on, 2022-11-14

[www.mabtech.com](http://www.mabtech.com)