

# Monoclonal antibody to Mouse Ig kappa chain

<b>Antibody:</b>	MT24
<b>Product code:</b>	3875-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Mouse immunoglobulins
<b>Isotype:</b>	Rat IgG1
<b>Specificity:</b>	Mouse Ig kappa light chain. No cross-reactivity (<0.1%) with Ig from rat, goat, rabbit, cow, sheep, horse or human.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	Antigen-specific ELISA and ELISpot
<b>Protocols:</b>	Please contact Mabtech for more information.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.