

Anti-human IgM mAb (MT22), ALP

Product code:	3880-9A-1000
Antibody:	MT22
Enzyme:	Highly purified enzyme from calf intestine with high activity especially suited for enzyme immunoassay purposes.
Conjugation:	The preparation has been thoroughly depleted of free ligand and free enzyme. The conjugate has high residual enzyme activity.
Quantity:	1 ml
Applications:	For quantification of IgM using e.g. ELISA. MAb MT22 is recommended as detection mAb in combination with coating mAbs MT11/12 (product code 3880-3).
Working dilution:	Supplied in 0.1 M Tris buffer with 0.002% Kathon CG. Dilute 1/1000 as a general guideline. Dilution factor may vary with differing assay conditions.
Storage:	Shipped at ambient temperature. On arrival store product at 4-8°C.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2022-11-22

www.mabtech.com