

Monoclonal Antibody to Mouse Ig lambda chain

BIOTINYLATED

| | |
|-----------------------|---|
| Antibody: | JC5-1 |
| Product code: | 3895-6-1000 |
| Quantity: | 1000 µg |
| Immunogen: | Mouse Ig lambda light chain |
| Isotype: | Rat IgG2b |
| Specificity: | Mouse Ig lambda light chain. No cross-reactivity (<0.1%) with Ig from rat, goat, rabbit, cow, sheep, horse or human. |
| Concentration: | Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide. |
| Purification: | Purified from <i>in vitro</i> cultures by protein G affinity chromatography. |
| Biotinylation: | Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin. |
| Storage: | Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use. |
| Applications: | Antigen-specific ELISA |
| Protocols: | Please contact Mabtech for more information. |

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com